


# BIOMVA

## Manual Vacuum Aspiration Kit



[www.portable-ultrasound-scanner.com](http://www.portable-ultrasound-scanner.com)

**Quality, Professional, Innovational**

# MVA kit

## Specifications

Name	Specification	Character	Application
Syringe	Length: 24.5cm Diameter: 3cm Capacity: 60ml Negative pressure value: 400mmhg	Easy operation Stable negative pressure Comfortable feeling Non-noisy	<ul style="list-style-type: none"> <li>• Abortion within 10 weeks</li> <li>• Clearing uterine surgery</li> <li>• Endometrial biopsy</li> </ul>
Cannula 4#	Length: 23cm Diameter: 4mm	End Yellow	
Cannula 5#	Length: 23cm Diameter: 5mm	End Green	
Cannula 6#	Length: 23cm Diameter: 6mm	End White	
Cannula 7#	Length: 23cm Diameter: 7mm	End White	



**Wuxi Biomedical Technology Co., Ltd.**      [www.portable-ultrasound-scanner.com](http://www.portable-ultrasound-scanner.com)

Add: Room 408, No.801 Hongqiao road, Wuxi, Jiangsu, China 214072

Tel: 86-510-85572789

Fax: 86-510-85572789

Email: [sales@cnbiomedical.com](mailto:sales@cnbiomedical.com)

Mobile phone: 86-13382885910



is a trade mark of Biomedical Technology Co., Ltd. Specifications subject to change without notice.